

## Product Datasheet

# Rabbit Anti-Sheep IgG (H&L) antibody (HRP) (orb1184977)

<b>Description</b>	Rabbit polyclonal antibody against Rabbit Anti-Sheep IgG (H&L) conjugated to HRP
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Sheep
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	ELISA, IHC, WB
<b>Immunogen</b>	Sheep IgG
<b>Preservatives</b>	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.05% Sodium Azide
<b>Concentration</b>	1mg/mL
<b>Storage</b>	Store at -20°C. Avoid freeze / thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Dilution Range</b>	E (1/5000 - 1/20000), WB (1/5000 - 1/20000), IH (1/500 - 1/1000), IC (1/500 - 1/1000)
<b>Expiration Date</b>	12 months from date of receipt.