

Product Datasheet

CD235 antibody (PE/Cy7) (orb1185092)

Description	Mouse monoclonal antibody against CD235 conjugated to PE/Cy7
Species/Host	Mouse
Reactivity	Human
Conjugation	PE/Cy7
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Storage	Store at 2~8°C. Do not freeze.
Note	For research use only
Isotype	Mouse IgG2b, κ
Uniprot ID	P06028 , P02724
Entrez	2993
Expiration Date	12 months from date of receipt.