

## Product Datasheet

### CD235 antibody (PE/Cy5) (orb1185093)

<b>Description</b>	Mouse monoclonal antibody against CD235 conjugated to PE/Cy5
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	PE/Cy5
<b>Tested Applications</b>	FC
<b>Preservatives</b>	PBS with 0.05% Proclin300, 1% rAlbumin
<b>Storage</b>	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
<b>Note</b>	For research use only
<b>Isotype</b>	Mouse IgG2b, κ
<b>Uniprot ID</b>	<a href="#">P06028</a> , <a href="#">P02724</a>
<b>Entrez</b>	<a href="#">2993</a>
<b>Expiration Date</b>	12 months from date of receipt.