

Product Datasheet

CD74 antibody (PE/Cy5.5) (orb1185119)

Description	Mouse monoclonal antibody against CD74 conjugated to PE/Cy5.5
Species/Host	Mouse
Reactivity	Human
Conjugation	PE/Cy5.5
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Storage	Store at 2~8°C. Do not freeze.
Note	For research use only
Isotype	Mouse IgG1, κ
Uniprot ID	P04233
Entrez	972
Expiration Date	12 months from date of receipt.