

Product Datasheet

CD183/CXCR3 antibody (Biotin) (orb1185194)

Description	Hamster monoclonal antibody against CD183/CXCR3 conjugated to Biotin
Species/Host	Hamster
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Concentration	0.5 mg/mL
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Armenian Hamster IgG
Uniprot ID	O88410
Entrez	12766
Expiration Date	12 months from date of receipt.