

Product Datasheet

CD16/32 antibody (PerCP/Cy5.5) (orb1185225)

Description	Rat monoclonal antibody against CD16/32 conjugated to PerCP/Cy5.5
Species/Host	Rat
Reactivity	Mouse
Conjugation	PerCP/Cy5.5
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Concentration	0.2 mg/mL
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Rat IgG2b, κ
Uniprot ID	P08508 , Q9ES92 , P08101 , Q7TMW9 , Q5D5J5 , Q5D5I8
Entrez	14131 , 14130
Expiration Date	12 months from date of receipt.