

Product Datasheet

CD122/IL-3RB antibody (APC) (orb1185314)

Description	Mouse monoclonal antibody against CD122/IL-3RB conjugated to APC
Species/Host	Mouse
Reactivity	Human
Conjugation	APC
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Mouse IgG1, κ
Uniprot ID	P14784
Entrez	3560
Expiration Date	12 months from date of receipt.