

Product Datasheet

CD32 antibody (APC) (orb1185355)

Description	Mouse monoclonal antibody against CD32 conjugated to APC
Species/Host	Mouse
Reactivity	Human
Conjugation	APC
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Storage	Store at 2~8°C. Do not freeze.
Note	For research use only
Isotype	Mouse IgG2b, κ
Uniprot ID	P31994
Entrez	2212
Expiration Date	12 months from date of receipt.