

Product Datasheet

CD28 antibody (PE/Cy7) (orb1185392)

Description	Hamster monoclonal antibody against CD28 conjugated to PE/Cy7
Species/Host	Hamster
Reactivity	Mouse
Conjugation	PE/Cy7
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Concentration	0.2 mg/mL
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Syrian Hamster IgG
Uniprot ID	P31041
Entrez	12487
Expiration Date	12 months from date of receipt.