

Product Datasheet

Ig light chain λ antibody (PE/Cy5.5) (orb1185460)

Description	Mouse monoclonal antibody against Ig light chain λ conjugated to PE/Cy5.5
Species/Host	Mouse
Reactivity	Human
Conjugation	PE/Cy5.5
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Mouse IgG2a, κ
Entrez	3535
Expiration Date	12 months from date of receipt.