

Product Datasheet

CD57 antibody (PE/Cy7) (orb1185529)

Description	Mouse monoclonal antibody against CD57 conjugated to PE/Cy7
Species/Host	Mouse
Reactivity	Human
Conjugation	PE/Cy7
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Mouse IgM
Uniprot ID	Q9P2W7
Entrez	27087
Expiration Date	12 months from date of receipt.