

## Product Datasheet

### CD16 antibody (PerCP/Cy5.5) (orb1185615)

<b>Description</b>	Mouse monoclonal antibody against CD16 conjugated to PerCP/Cy5.5
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	PerCP/Cy5.5
<b>Tested Applications</b>	FC
<b>Preservatives</b>	PBS with 0.05% Proclin300, 1% rAlbumin
<b>Storage</b>	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
<b>Note</b>	For research use only
<b>Isotype</b>	Mouse IgG1, κ
<b>Uniprot ID</b>	<a href="#">P08637</a>
<b>Entrez</b>	<a href="#">2214</a>
<b>Expiration Date</b>	12 months from date of receipt.