

## Product Datasheet

### Ly-6G/Ly-6C (Gr-1) antibody (PE/Cy5.5) (orb1185687)

<b>Description</b>	Rat monoclonal antibody against Ly-6G/Ly-6C (Gr-1) conjugated to PE/Cy5.5
<b>Species/Host</b>	Rat
<b>Reactivity</b>	Mouse
<b>Conjugation</b>	PE/Cy5.5
<b>Tested Applications</b>	FC
<b>Preservatives</b>	PBS with 0.05% Proclin300, 1% rAlbumin
<b>Concentration</b>	0.2 mg/mL
<b>Storage</b>	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
<b>Note</b>	For research use only
<b>Isotype</b>	Rat IgG2b, $\kappa$
<b>Uniprot ID</b>	<a href="#">POCW03</a> , <a href="#">P35461</a> , <a href="#">POCW02</a>
<b>Entrez</b>	<a href="#">17067</a> , <a href="#">546644</a>
<b>Expiration Date</b>	12 months from date of receipt.