

Product Datasheet

Ly6G antibody (PE/Cy5.5) (orb1185703)

Description	Rat monoclonal antibody against Ly6G conjugated to PE/Cy5.5
Species/Host	Rat
Reactivity	Mouse
Conjugation	PE/Cy5.5
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Concentration	0.2 mg/mL
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Rat IgG2a, κ
Uniprot ID	P35461
Entrez	546644
Expiration Date	12 months from date of receipt.