

Product Datasheet

CD273/PD-L2 antibody (PE) (orb1185860)

Description	Mouse monoclonal antibody against CD273/PD-L2 conjugated to PE
Species/Host	Mouse
Reactivity	Human
Conjugation	PE
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Mouse IgG2a, κ
Uniprot ID	Q9BQ51
Entrez	80380
Expiration Date	12 months from date of receipt.