

Product Datasheet

CD86 antibody (PE/Cy7) (orb1186139)

Description	Rat monoclonal antibody against CD86 conjugated to PE/Cy7
Species/Host	Rat
Reactivity	Mouse
Conjugation	PE/Cy7
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Concentration	0.2 mg/mL
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Rat IgG2a, κ
Uniprot ID	P42082
Entrez	12524
Expiration Date	12 months from date of receipt.