

Product Datasheet

Ly-6G/Ly-6C (Gr-1) antibody (PE/Cy5) (orb1186236)

Description	Rat monoclonal antibody against Ly-6G/Ly-6C (Gr-1) conjugated to PE/Cy5
Species/Host	Rat
Reactivity	Mouse
Conjugation	PE/Cy5
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Concentration	0.2 mg/mL
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Rat IgG2b, κ
Uniprot ID	POCW03 , P35461 , POCW02
Entrez	17067 , 546644
Expiration Date	12 months from date of receipt.