

Product Datasheet

CD11a antibody (Biotin) (orb1186270)

Description	Rat monoclonal antibody against CD11a conjugated to Biotin
Species/Host	Rat
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Concentration	0.5 mg/mL
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Rat IgG2b, κ
Uniprot ID	P24063
Entrez	16408
Expiration Date	12 months from date of receipt.