

Product Datasheet

CD122 antibody (APC) (orb1186273)

Description	Rat monoclonal antibody against CD122 conjugated to APC
Species/Host	Rat
Reactivity	Mouse
Conjugation	APC
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Concentration	0.2 mg/mL
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Rat IgG2a, κ
Uniprot ID	P16297
Entrez	16185
Expiration Date	12 months from date of receipt.