

Product Datasheet

CD86 antibody (PE/Cy7) (orb1186321)

| | |
|----------------------------|---|
| Description | Rat monoclonal antibody against CD86 conjugated to PE/Cy7 |
| Species/Host | Rat |
| Reactivity | Mouse |
| Conjugation | PE/Cy7 |
| Tested Applications | FC |
| Preservatives | PBS with 0.05% Proclin300, 1% rAlbumin |
| Storage | Store at 2~8°C and protected from prolonged exposure to light. Do not freeze. |
| Note | For research use only |
| Isotype | Rat IgG2a, κ |
| Uniprot ID | P42082 |
| Entrez | 12524 |
| Expiration Date | 12 months from date of receipt. |