

Product Datasheet

CD71 antibody (Biotin) (orb1186453)

Description	Rat monoclonal antibody against CD71 conjugated to Biotin
Species/Host	Rat
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Concentration	0.5 mg/mL
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Rat IgG2a, κ
Uniprot ID	Q62351
Entrez	22042
Expiration Date	12 months from date of receipt.