

Product Datasheet

Rat IgG2b, κ Isotype Control (PE/Cy5.5) (orb1186474)

| | |
|------------------------|---|
| Description | Rat IgG2b, κ Isotype Control (PE/Cy5.5) |
| Species/Host | Rat |
| Conjugation | PE/Cy5.5 |
| Preservatives | PBS with 0.05% Proclin300, 1% rAlbumin |
| Concentration | 0.5 mg/mL |
| Storage | Store at 2~8°C and protected from prolonged exposure to light. Do not freeze. |
| Note | For research use only |
| Isotype | IgG2b, κ |
| Expiration Date | 12 months from date of receipt. |