

Product Datasheet

Rat IgG2a, κ Isotype Control (APC/Cy5.5) (orb1186481)

Description	Rat IgG2a, κ Isotype Control (APC/Cy5.5)
Species/Host	Rat
Conjugation	APC/Cy5.5
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Concentration	0.5 mg/mL
Storage	Store at 2~8°C. Do not freeze.
Note	For research use only
Isotype	IgG2a, κ
Expiration Date	12 months from date of receipt.