

Product Datasheet

Rat IgG1, κ Isotype Control (PE/Cy5.5) (orb1186492)

Description	Rat IgG1, κ Isotype Control (PE/Cy5.5)
Species/Host	Rat
Conjugation	PE/Cy5.5
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Concentration	0.5 mg/mL
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	IgG1, κ
Expiration Date	12 months from date of receipt.