

Product Datasheet

Mouse IgG2b, κ Isotype Control (PE/Cy7) (orb1186502)

Description	Mouse IgG2b, κ Isotype Control (PE/Cy7)
Species/Host	Mouse
Conjugation	PE/Cy7
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Concentration	0.5 mg/mL
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	IgG2b, κ
Expiration Date	12 months from date of receipt.