

## Product Datasheet

### Mouse IgG2b, $\kappa$ Isotype Control (PE/Cy5) (orb1186503)

<b>Description</b>	Mouse IgG2b, $\kappa$ Isotype Control (PE/Cy5)
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	PE/Cy5
<b>Preservatives</b>	PBS with 0.05% Proclin300, 1% rAlbumin
<b>Concentration</b>	0.5 mg/mL
<b>Storage</b>	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2b, $\kappa$
<b>Expiration Date</b>	12 months from date of receipt.