

Product Datasheet

Mouse IgG2a, κ Isotype Control (PerCP/Cy5.5) (orb1186509)

Description	Mouse IgG2a, κ Isotype Control (PerCP/Cy5.5)
Species/Host	Mouse
Conjugation	PerCP/Cy5.5
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Concentration	0.5 mg/mL
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	IgG2a, κ
Expiration Date	12 months from date of receipt.