

Product Datasheet

Rat IgM, κ Isotype Control (PerCP/Cy5.5) (orb1186536)

Description	Rat IgM, κ Isotype Control (PerCP/Cy5.5)
Species/Host	Rat
Conjugation	PerCP/Cy5.5
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	IgM, κ
Expiration Date	12 months from date of receipt.