

## Product Datasheet

### Rat IgG2a, $\kappa$ Isotype Control (PerCP/Cy5.5) (orb1186579)

|                        |                                                   |
|------------------------|---------------------------------------------------|
| <b>Description</b>     | Rat IgG2a, $\kappa$ Isotype Control (PerCP/Cy5.5) |
| <b>Species/Host</b>    | Rat                                               |
| <b>Conjugation</b>     | PerCP/Cy5.5                                       |
| <b>Preservatives</b>   | PBS with 0.05% Proclin300,<br>1% rAlbumin         |
| <b>Concentration</b>   | 0.2 mg/mL                                         |
| <b>Storage</b>         | Store at 2~8°C. Do not freeze.                    |
| <b>Note</b>            | For research use only                             |
| <b>Isotype</b>         | IgG2a, $\kappa$                                   |
| <b>Expiration Date</b> | 12 months from date of receipt.                   |