

## Product Datasheet

# Mouse IgG2b, $\kappa$ Isotype Control (PE/Cy5.5) (orb1186598)

|                        |   |
|------------------------|---|
| <b>Description</b>     | Mouse IgG2b, $\kappa$ Isotype Control (PE/Cy5.5)                              |
| <b>Species/Host</b>    | Mouse   |
| <b>Conjugation</b>     | PE/Cy5.5  |
| <b>Preservatives</b>   | PBS with 0.05% Proclin300, 1% rAlbumin  |
| <b>Concentration</b>   | 0.2 mg/mL   |
| <b>Storage</b>         | Store at 2~8°C and protected from prolonged exposure to light. Do not freeze. |
| <b>Note</b>            | For research use only   |
| <b>Isotype</b>         | IgG2b, $\kappa$   |
| <b>Expiration Date</b> | 12 months from date of receipt.   |