

## Product Datasheet

# Rat IgM, $\kappa$ Isotype Control (PerCP/Cy5.5) (orb1186633)

<b>Description</b>	Rat IgM, $\kappa$ Isotype Control (PerCP/Cy5.5)
<b>Species/Host</b>	Rat
<b>Conjugation</b>	PerCP/Cy5.5
<b>Tested Applications</b>	FC
<b>Preservatives</b>	PBS with 0.05% Proclin300, 1% rAlbumin
<b>Concentration</b>	0.2 mg/mL
<b>Storage</b>	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
<b>Note</b>	For research use only
<b>Isotype</b>	IgM, $\kappa$
<b>Expiration Date</b>	12 months from date of receipt.