

Product Datasheet

Mouse IgM, κ Isotype Control (PE/Cy5.5) (orb1186643)

Description	Mouse IgM, κ Isotype Control (PE/Cy5.5)
Species/Host	Mouse
Conjugation	PE/Cy5.5
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Concentration	0.2 mg/mL
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	IgM, κ
Expiration Date	12 months from date of receipt.