

Product Datasheet

CD15 antibody (Biotin) (orb1186778)

Description	Mouse monoclonal antibody against CD15 conjugated to Biotin
Species/Host	Mouse
Reactivity	Human
Conjugation	Biotin
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Concentration	0.5 mg/mL
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Mouse IgM
Uniprot ID	P22083
Entrez	2526
Expiration Date	12 months from date of receipt.