

Product Datasheet

CD11a antibody (Biotin) (orb1186978)

| | |
|----------------------------|---|
| Description | Mouse monoclonal antibody against CD11a conjugated to Biotin |
| Species/Host | Mouse |
| Reactivity | Human |
| Conjugation | Biotin |
| Tested Applications | FC |
| Preservatives | PBS with 0.05% Proclin300, 1% rAlbumin |
| Concentration | 0.5 mg/mL |
| Storage | Store at 2~8°C and protected from prolonged exposure to light. Do not freeze. |
| Note | For research use only |
| Isotype | Mouse IgG1, κ |
| Uniprot ID | P20701 |
| Entrez | 3683 |
| Expiration Date | 12 months from date of receipt. |