

Product Datasheet

F4/80 antibody (PerCP/Cy5.5) (orb1187425)

Description	Rat monoclonal antibody against F4/80 conjugated to PerCP/Cy5.5
Species/Host	Rat
Reactivity	Mouse
Conjugation	PerCP/Cy5.5
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Rat IgG2b, κ
Uniprot ID	Q61549
Entrez	13733
Expiration Date	12 months from date of receipt.