

Product Datasheet

Rat IgM, κ Isotype Control (APC) (orb1187486)

Description	Rat IgM, κ Isotype Control (APC)
Species/Host	Rat
Conjugation	APC
Tested Applications	FC
Preservatives	PBS with 0.05% Proclin300, 1% rAlbumin
Concentration	0.2 mg/mL
Storage	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	IgM, κ
Expiration Date	12 months from date of receipt.