

## Product Datasheet

### Rat IgG2b, $\kappa$ Isotype Control (orb1187535)

<b>Description</b>	Rat IgG2b, $\kappa$ Isotype Control
<b>Species/Host</b>	Rat
<b>Conjugation</b>	Unconjugated
<b>Preservatives</b>	PBS with 0.05% Proclin300
<b>Concentration</b>	0.2 mg/mL
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2b, $\kappa$
<b>Expiration Date</b>	12 months from date of receipt.