

Product Datasheet

Mouse IgG2b, κ Isotype Control (orb1187551)

Description	Mouse IgG2b, κ Isotype Control
Species/Host	Mouse
Conjugation	Unconjugated
Preservatives	0.2 μ m filtered in PBS, pH7.2. Contains no stabilizers or preservatives. Endotoxin level is 2 EU/mg (0.002EU/ μ g) as Determined by LAL gel clotting assay. Purity is > 95% as determined by SDS-PAGE.
Concentration	0.2 mg/mL
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG2b, κ
Expiration Date	12 months from date of receipt.