

Product Datasheet

Mouse IgG1, κ Isotype Control (orb1187560)

Description	Mouse IgG1, κ Isotype Control
Species/Host	Mouse
Conjugation	Unconjugated
Preservatives	PBS with 0.05% Proclin300
Concentration	0.2 mg/mL
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1, κ
Expiration Date	12 months from date of receipt.