

Product Datasheet

Rat IgM, κ Isotype Control (Biotin) (orb1187569)

| | |
|----------------------------|---|
| Description | Rat IgM, κ Isotype Control (Biotin) |
| Species/Host | Rat |
| Conjugation | Biotin |
| Tested Applications | FC |
| Preservatives | PBS with 0.05% Proclin300, 1% rAlbumin |
| Concentration | 0.5 mg/mL |
| Storage | Store at 2~8°C and protected from prolonged exposure to light. Do not freeze. |
| Note | For research use only |
| Isotype | IgM, κ |
| Expiration Date | 12 months from date of receipt. |