

Product Datasheet

HA-tag antibody (orb1187588)

Description	Mouse monoclonal antibody to HA-tag
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	IF, IP, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence of HA-tag. The exact sequence is proprietary.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/mL
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Dilution Range	WB (1/1000 - 1/3000), IF/IC (1/100 - 1/200), IP (1/50 - 1/200)
Expiration Date	12 months from date of receipt.