

## Product Datasheet

# Hepatitis B Virus Surface Antigen antibody (Biotin) (orb1187600)

<b>Description</b>	Mouse polyclonal antibody against Hepatitis B Virus Surface Antigen conjugated to Biotin
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Biotin
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Recombinant Hepatitis B Virus Surface Antigen
<b>Preservatives</b>	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.
<b>Concentration</b>	1mg/mL
<b>Storage</b>	Store at -20°C. Avoid repeated freeze / thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Dilution Range</b>	ELISA: Use at an assay dependent dilution.
<b>Expiration Date</b>	12 months from date of receipt.