

Product Datasheet

Hepatitis B Virus Core Antigen antibody (FITC) (orb1187609)

| | |
|----------------------------|--|
| Description | Mouse monoclonal antibody against Hepatitis B Virus Core Antigen conjugated to FITC |
| Species/Host | Mouse |
| Reactivity | Virus |
| Conjugation | FITC |
| Tested Applications | ELISA |
| Immunogen | Recombinant Hepatitis B Virus Core Antigen |
| Preservatives | Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide. |
| Concentration | 1mg/mL |
| Storage | Store at -20°C. Avoid repeated freeze / thaw cycles. |
| Note | For research use only |
| Clonality | Monoclonal |
| Dilution Range | ELISA: Use at an assay dependent dilution. |
| Expiration Date | 12 months from date of receipt. |