

Product Datasheet

Hepatitis B Virus Core Antigen antibody (FITC) (orb1187610)

Description	Goat monoclonal antibody against Hepatitis B Virus Core Antigen conjugated to FITC
Species/Host	Goat
Reactivity	Virus
Conjugation	FITC
Tested Applications	ELISA
Immunogen	Recombinant Hepatitis B Virus Core Antigen
Preservatives	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.
Concentration	1mg/mL
Storage	Store at -20°C. Avoid repeated freeze / thaw cycles.
Note	For research use only
Clonality	Monoclonal
Dilution Range	ELISA: Use at an assay dependent dilution.
Expiration Date	12 months from date of receipt.