

Product Datasheet

Nano-tag antibody (orb1187713)

Description	Mouse monoclonal antibody to Nano-tag
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence of Nano-tag. The exact sequence is proprietary.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/mL
Storage	Store at -20°C. Avoid repeated freeze / thaw cycles.
Note	For research use only
Clonality	Monoclonal
Dilution Range	WB (1/1000 - 1/2000), IH (1/200 - 1/500)
Expiration Date	12 months from date of receipt.