

Product Datasheet

His-tag antibody (orb1187751)

Description	Mouse monoclonal antibody to His-tag
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH-conjugated synthetic peptide HHHHHH.
Preservatives	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.
Concentration	1mg/mL
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Dilution Range	ELISA: Use at an assay dependent dilution.
Expiration Date	12 months from date of receipt.