

## Product Datasheet

### SRT-tag antibody (orb1187800)

|                            |                                                                                                           |
|----------------------------|-----------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to SRT-tag                                                                      |
| <b>Species/Host</b>        | Mouse                                                                                                     |
| <b>Conjugation</b>         | Unconjugated                                                                                              |
| <b>Tested Applications</b> | WB                                                                                                        |
| <b>Immunogen</b>           | Recombinant protein corresponding to SRT-tag.                                                             |
| <b>Preservatives</b>       | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide. |
| <b>Concentration</b>       | 1mg/mL                                                                                                    |
| <b>Storage</b>             | Store at -20°C. Avoid repeated freeze / thaw cycles.                                                      |
| <b>Note</b>                | For research use only                                                                                     |
| <b>Clonality</b>           | Monoclonal                                                                                                |
| <b>Dilution Range</b>      | WB (1/2000 - 1/5000)                                                                                      |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                           |