

Product Datasheet

HSV-tag antibody (orb1187810)

| | |
|----------------------------|---|
| Description | Mouse monoclonal antibody to HSV-tag |
| Species/Host | Mouse |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | KLH-conjugated synthetic peptide encompassing a sequence of HSV-tag. The exact sequence is proprietary. |
| Preservatives | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide. |
| Concentration | 1mg/mL |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Clonality | Monoclonal |
| Dilution Range | WB (1/5000 - 1/10000) |
| Expiration Date | 12 months from date of receipt. |