

Product Datasheet

FLAG-tag antibody (HRP) (orb1187817)

| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody against FLAG-tag conjugated to HRP |
| Species/Host | Rabbit |
| Conjugation | HRP |
| Tested Applications | WB |
| Immunogen | KLH-conjugated synthetic peptide encompassing a sequence of FLAG-tag. The exact sequence is proprietary. |
| Preservatives | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide. |
| Concentration | 1mg/mL |
| Storage | Store at -20°C. Avoid repeated freeze / thaw cycles. |
| Note | For research use only |
| Clonality | Polyclonal |
| Dilution Range | WB (1/5000 - 1/10000) |
| Expiration Date | 12 months from date of receipt. |