

## Product Datasheet

### GST-tag antibody (orb1187826)

<b>Description</b>	Rabbit polyclonal antibody to GST-tag
<b>Species/Host</b>	Rabbit
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IP, WB
<b>Immunogen</b>	KLH-conjugated synthetic peptide SVDLMKKRLENGDD.
<b>Preservatives</b>	Liquid in PBS, pH 7.3, 0.2% rAlbumin, and 0.02% sodium azide.
<b>Concentration</b>	1mg/mL
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Dilution Range</b>	WB (1/500 - 1/1000), IP (1/10 - 1/100)
<b>Expiration Date</b>	12 months from date of receipt.